


OLYMPUS

Your Vision, Our Future

STONE MANAGEMENT

UROPASS – URETERAL ACCESS SHEATH

Access with Confidence.




UROPASS – THE SOLUTION FOR URETERAL ACCESS

Access Should Never Be the Challenge

Protecting the patient's anatomy is not just for the most challenging of stone cases. Accessing the kidney and multiple scope passes can become a challenge even in what seems like a straightforward case. By utilizing the UroPass Access Sheath, experience tells us that accelerated results can be achieved while maintaining the integrity of the anatomy.

- Complete portfolio: Availability of different diameters from 10/12 to 13/15 Fr.
- Protection against elevated renal pressures during flexible ureteroscopy.
- Reduced procedure times by enabling multiple passes while protecting the ureter.

Patented clip and funnel design ensures sheath dilator tip position is maintained during placement




Hydrophilic coating reduces friction and facilitates passage

Olympus Quality for Every Situation

The UroPass Access Sheath platform is designed for compatibility with our scopes, leveraging years of expertise in access and visualization. You can think of our UroPass Access Sheath as a scope within a scope, providing reach, safety, confidence, and design to help you deliver better outcomes.


Radiopaque tip and sheath provide easy viewing of the placement location



Stainless-steel-reinforced sheath prevents kinking to avoid the ureteroscope getting stuck



Atraumatic tip dilation from 6 to 12 Fr



14189

287

14190

ORDERING INFORMATION

Ureteral Access Sheaths			
Cat. Number	Diameter	Length	Pieces/Box
EG61024BX	10/12 Fr	24 cm	5
EG61038BX	10/12 Fr	38 cm	5
EG61046BX	10/12 Fr	46 cm	5
EG61054BX	10/12 Fr	54 cm	5
EG61124BX	11/13 Fr	24 cm	5
EG61138BX	11/13 Fr	38 cm	5
EG61146BX	11/13 Fr	46 cm	5
EG61154BX	11/13 Fr	54 cm	5

Ureteral Access Sheaths, Continued			
Cat. Number	Diameter	Length	Pieces/Box
EG61224BX	12/14 Fr	24 cm	5
EG61238BX	12/14 Fr	38 cm	5
EG61246BX	12/14 Fr	46 cm	5
EG61254BX	12/14 Fr	54 cm	5
EG61324BX	13/15 Fr	24 cm	5
EG61338BX	13/15 Fr	38 cm	5
EG61346BX	13/15 Fr	46 cm	5
EG61354BX	13/15 Fr	54 cm	5

Olympus Stone Management Solutions

URF-V2 – superslim 8.4 Fr uretero-reno videoscope



8857

Ultra-Catch No-Tip Nitinol Basket



10961

Ureteral Stents



14195

URF-P6 – superslim 7.95 Fr uretero-reno fiberscope



4169

Irri-Flo II – irrigation pump for improved flow and visualization



5753

Guidewires



10965

Specifications, design, and accessories are subject to change without any notice or obligation on the part of the manufacturer.

OLYMPUS

OLYMPUS EUROPA SE & CO. KG
 Postbox 10 49 08, 20034 Hamburg, Germany
 Wendenstrasse 14-18, 20097 Hamburg, Germany
 Phone: +49 40 23773-0, Fax: +49 40 233765
 www.olympus-europa.com



E0492343 · 1.500 · 06/15 · PR · HB