

## PK Morcellator

Powering Gynecology



### Peeling Tip

Integrated peeling tip promotes peeling on target tissue and avoids tissue coring throughout the procedure.

### Reduced Surgical Smoke\*

The smoke management of the PK Morcellator maintains clear visibility of the morcellator, the target tissue, and collateral tissue.

### Faster Morcellation Speed\*

PK Morcellator completes morcellation at a fast rate.



\* Compared to PKS PlasmaSORD

# PK Morcellator

## Ordering Information

Order Number	Product Description	Diameter	Length	Pieces
<b>PK TECHNOLOGY Instruments</b>				
<b>EGPK-CF0533</b>	PK Cutting Forceps	5 mm	33 cm	5/box
<b>WA90300A</b>	PK Lap Loop	5 mm	33 cm	5/box
<b>WA90200A</b>	PK Morcellator	15 mm	17 cm	5/box
<b>EGPK-NE0533</b>	PK Needle	5 mm	33 cm	5/box
<b>EGPK-JH0533</b>	PK J-Hook	5 mm	33 cm	5/box
<b>EGPK-SP0533</b>	PK Spatula	5 mm	33 cm	5/box

## Generators and Accessories

<b>WB91051W</b>	ESG-400 electrosurgical generator
<b>WA95623A</b>	Power cable UK plug, 4.5 m (1x)
<b>WA95621A</b>	Power cable Europlug, 4.5 m (1x)

## Optional Equipment

<b>N3635730</b>	Energy cart, TC-E400
<b>WB50403W</b>	Foot switch, single (bipolar) pedal
<b>N3809530</b>	Adaptor for the connection to UHI-3 or UHI-4 insufflator
<b>E0421213</b>	Neutral electrode cable (reusable)

Specifications, design, and accessories are subject to change without any notice or obligation on the part of the manufacturer.

**OLYMPUS**

**OLYMPUS EUROPA SE & CO. KG**  
Postbox 10 49 08, 20034 Hamburg, Germany  
Wendenstrasse 14-18, 20097 Hamburg, Germany  
Phone: +49 40 23773-0, Fax: +49 40 233765  
www.olympus-europa.com