


PK Lap Loop

Powering Gynecology



Patient Safety


- Blue-coated wire for optimal visibility
- White ceramic tip helps to identify if bowel is caught in loop

Efficiency

- Fast, clean, and immediate bipolar cutting during LSH procedures¹

Easy and Convenient Handling and Placement

- Large loop for big uterus
- Pop-up diamond-shaped loop opens up automatically in the abdomen
- Consistent loop shape
- Ergonomic handle design



¹Erian J, et al.: "One hundred cases of LSH using the PK and Lap Loop systems"; J Min Invasive Gyn 2005;12:365-369; Competitor: Medsys.

PK Lap Loop

Ordering Information

Order Number	Product Description	Diameter	Length	Pieces
PK TECHNOLOGY Instruments				
EGPK-CF0533	PK Cutting Forceps	5 mm	33 cm	5/box
WA90300A	PK Lap Loop	5 mm	33 cm	5/box
WA90200A	PK Morcellator	15 mm	17 cm	5/box
EGPK-NE0533	PK Needle	5 mm	33 cm	5/box
EGPK-JH0533	PK J-Hook	5 mm	33 cm	5/box
EGPK-SP0533	PK Spatula	5 mm	33 cm	5/box

Generators and Accessories

WB91051W	ESG-400 electrosurgical generator
WA95623A	Power cable UK plug, 4.5 m (1x)
WA95621A	Power cable Europlug, 4.5 m (1x)

Optional Equipment

N3635730	Energy cart, TC-E400
WB50403W	Foot switch, single (bipolar) pedal
N3809530	Adaptor for the connection to UHI-3 or UHI-4 insufflator
E0421213	Neutral electrode cable (reusable)

Specifications, design, and accessories are subject to change without any notice or obligation on the part of the manufacturer.

OLYMPUS

OLYMPUS EUROPA SE & CO. KG

Postbox 10 49 08, 20034 Hamburg, Germany
Wendenstrasse 14-18, 20097 Hamburg, Germany
Phone: +49 40 23773-0, Fax: +49 40 233765
www.olympus-europa.com